

Airnet® II

2-Channel Particle Sensor

Easy to install, compact particle counters meeting the latest international regulations

Without measurement there is no control

The Airnet® II particle sensor makes it easy and cost-effective to monitor your cleanroom. This compact particle sensor offers unparalleled performance and data transmission capabilities while meeting ISO 21501-4 and ISO 14644-1:2015.

With a simple installation and versatile power options, the unit is fully configurable to accept distributed power from an in-house system, local power plug-in, or Power over Ethernet (PoE). Communication capabilities include Ethernet to interface with Pharmaceutical Net Pro, Facility Net, or FacilityPro® software, Modbus communications, or optional 4-20 mA output.

Data integrity is maintained through the use of a data queue feature that continues to gather data even if network communication is lost.

To ensure proper flow conditions and vacuum system operation, the Airnet II sensors incorporate a dynamic flow sensing system that alarms with a 15% change in flow conditions.

For applications where decontamination using Vaporized Hydrogen Peroxide (VHP) is required, an optional VHP-compatible unit is available for a simple installation without complex valving (Airnet 510XR).

AVAILABLE WITH
ISO 17025
CALIBRATION
for select parameters



BENEFITS

- Proven technology provides reliable and accurate data
- Immediate reaction to particle contamination events
- Low-cost solution for multipoint monitoring
- Interfaces with Facility Net, Pharmaceutical Net Pro and FacilityPro Software for comprehensive management of cleanroom conditions
- Small size and flexible mounting options make it easy to install in cleanrooms and mini-environments
- Laser diode drastically reduces the need for maintenance and extends product lifecycle
- Automatic laser shutdown reduces laser failures
- Data queue maintains data integrity when communication is lost
- Optional 4-20 mA output for integration with existing systems
- Optional XR coating protects sensors against corrosive or oxidizing vapors in VHP sterilization processes

FEATURES

- 2 channels
- 0.5 and 5.0 μm size range, suitable for pharmaceutical applications
- 1.0 CFM flowrate
- Interfaces with Modbus communications and optional 4-20 mA output
- Chemical-resistant polycarbonate enclosure
- Low sample point cost
- Small enough for use in remote locations
- ISO 14644-1:2015 compliant

APPLICATIONS

- Cleanroom monitoring
- Dedicated monitoring of critical locations
- Trend analysis
- Statistical process control
- Multi-location monitoring
- Isolator monitoring

Airnet® II

2 Channel Particle Sensor

Specifications

	301	501/501A	510/510XR
Size Range	0.3, 0.5 µm	0.5, 5.0 µm ¹	0.5, 5.0 µm
Maximum Concentration ²	5,695,168/ft ³	9,578,238/ft ³	957,824/ft ³
Flow Rate	0.1 CFM (2.8 LPM)		1.0 CFM (28.3 LPM)
Zero Count	≤ 70.7 counts/m ³		≤ 7.07 counts/m ³
Vacuum Source	> 11 in Hg (> 375 mBar) below atm pressure		> 15 in Hg (> 510 mBar) below atm pressure
Counting Efficiency	50% ± 20% for most-sensitive channel. Meets ISO 21501-4. 100% ± 10% at 1.5 to 2.0 times channel one size. Meets ISO 21501-4.		
Laser Source	Diode		
Laser Classification	Class 1 per EN 60825 (Internally, a Class IIIB laser is used, per EN 60825)		
Exterior Surface	Polycarbonate		
Dimensions (H x W x L)	5.3 x 3.6 x 3.8 in (13.5 x 8.9 x 9.6 cm)		
Weight	1.6 lb (0.73 kg)		
Sample Probe or Tubing	1/4-in ID		
Flow System	External vacuum 1/4-in connection; automatic laser shutoff and alarm on 15% flow variation		
Power	24 VDC (0.5 A) 100 – 240 VAC ± 10%, 50 – 60 Hz, Power Supply (optional), or use of Power over Ethernet 48 VDC via a PoE router		
Communication Connectors	Ethernet (Particle Measuring Systems proprietary protocol, Modbus TCP) RS-232 (configuration and diagnostic tool only, no data) 4-20 mA, optional (3 output channels: 2 particle data, 1 instrument status)		
Status Indicators	Programmable status (two-color LED), Activity (one-color LED)		
Calibration	Calibration materials used are traceable to the National Institute of Standards and Technology (NIST) and meet ISO 21501-4 requirements.		
Environment	Temperature: 4 – 35 °C (39 – 95 °F), 5 – 95%; non-condensing relative humidity		
Complies with	EU RoHS, ISO 21501-4, ISO 14644-1:2015		

¹ Airnet II 501A has 0.5 and 1.0 µm channels.

² 10% coincidence loss at maximum concentration.

HEADQUARTERS

5475 Airport Blvd
Boulder, Colorado 80301 USA
T: +1 303 443 7100, +1 800 238 1801

Instrument Service & Support
T: +1 800 557 6363

Customer Response Center
T: +1 877 475 3317
E: info@pmeasuring.com

www.pmeasuring.com
info@pmeasuring.com



**PARTICLE
MEASURING
SYSTEMS®**
a spectris company

GLOBAL OFFICES

AUSTRIA
T: +43 512 390 500
E: pmsaustria@pmeasuring.com

BENELUX
T: +32 10 23 71 56
E: pmsbelgium@pmeasuring.com

BRAZIL
T: +55 11 5188 8227
E: pmsbrazil@pmeasuring.com

CHINA
T: +86 21 6113 3600
E: pmschina@pmeasuring.com

FRANCE
T: +33(0)1 60 10 32 96
E: pmsfrance@pmeasuring.com

GERMANY
T: +49 6151 6671 632
E: pmsgermany@pmeasuring.com

ITALY
T: +39 06 9053 0130
E: pmsrl@pmeasuring.com

JAPAN
T: +81 44 589 3498
E: pmsjapan@pmeasuring.com

KOREA
T: +82 31 286 5790
E: pmskorea@pmeasuring.com

MEXICO
T: +52 55 2271 5106
E: pmsmexico@pmeasuring.com

NORDIC
T: +45 707 028 55
E: pmsnordic@pmeasuring.com

PUERTO RICO
T: +1 787 718 9096
E: pmspuertorico@pmeasuring.com

SINGAPORE
T: +65 6496 0330
E: pmsingapore@pmeasuring.com

SWITZERLAND
T: +41 71 987 01 01
E: pmsswitzerland@pmeasuring.com

TAIWAN
T: 886-3-5525300 Ext: 301
E: pms taiwan@pmeasuring.com